

# Amber®

*PushClick Technology*



## *Technical Specifications*

<b>Device Type</b>	Disposable Autoinjector
<b>Primary Container</b>	1mL Long Pre-filled Syringe   FNS or RNS
<b>Dosage Type</b>	Fixed   Full Contents
<b>Needle Extension</b>	6 - 8 mm   Subcutaneous
<b>Needle Insertion</b>	Automatic
<b>Needle Removal</b>	Manual
<b>Activation Method</b>	Needle Cover

# Handling Instructions



Remove cap



To activate, press & hold against injection site



Remove device from injection site & discard in sharps container

## Key Features

- ✓ Ergonomic industrial design
- ✓ Permanently hidden needle
- ✓ Needle cover activation
- ✓ Audible, visual & tactile feedback