



Amber®

PushClick Technology

Volume 

Viscosity 

Device Type

Disposable Auto Injector

Primary Container

1mL Long Pre-filled Syringe | FNS or RNS

Dosage Type

Fixed | Full Contents

Route of Administration

Subcutaneous

Needle Insertion

Automatic

Activation Method

Needle Cover

Injection Feedback

Audible click at beginning & end of injection, viewing window

Handling Instructions



1
Remove cap



2
Press & hold against injection site
to activate



+
Discard in sharps container

Key Features

- ✓ Simple 2-step operation
- ✓ Needle cover activation
- ✓ Automatic needle insertion
- ✓ Audible, visual & tactile feedback
- ✓ Permanently hidden needle
- ✓ One-handed operation
- ✓ Pull- or twist-off cap
- ✓ Ergonomic industrial design
- ✓ Passive locking needle cover