



Bertha®

For high volume & viscosity

Volume Viscosity **Device Type**

Disposable Auto Injector

Primary Container

1mL and 2.25mL Syringe with Staked Needle

Dosage Type

Fixed Dose

Route of Administration

Subcutaneous

Needle Insertion

Manual

Activation Method

Needle Cover

Injection Feedback

Continuous clicks during injection, visual feedback ribbon & viewing window

Handling Instructions



Remove cap



Press & hold against injection site to activate



Discard in sharps container

Key Features

- ✓ Simple 2-step operation
- ✓ Needle cover activation
- ✓ Manual needle insertion
- ✓ Continuous audible clicks & visual feedback
- ✓ One-handed operation
- ✓ Permanently hidden needle
- ✓ Passive locking needle cover
- ✓ Preconfigured or customized to client specifications
- ✓ Compatible with Rigid Needle Shield