



Bertha®

For high volume & viscosity

Volume 

Viscosity 

Device Type	Disposable Auto Injector
Primary Container	1mL or 2.25mL Syringe with Staked Needle
Dosage Type	Fixed Dose
Route of Administration	Subcutaneous
Needle Insertion	Manual
Activation Method	Needle Cover
Injection Feedback	Continuous clicks during injection, visual feedback ribbon & viewing window

Handling Instructions



1
Remove cap



2
Press & hold against injection site to activate



+Discard in sharps container

Key Features

- ✓ Simple 2-step operation
- ✓ Needle cover activation
- ✓ Manual needle insertion
- ✓ Continuous audible clicks & visual feedback
- ✓ One-handed operation
- ✓ Permanently hidden needle
- ✓ Passive locking needle cover
- ✓ Preconfigured or customized to client specifications
- ✓ Compatible with Rigid Needle Shield