



DAI[®] 2

Revolutionary 2-Step Solution

Volume 

Viscosity 

Device Type

Disposable Auto Injector

Primary Container

1mL Long Syringe with Staked Needle

Dosage Type

Fixed | Full Contents

Route of Administration

Subcutaneous

Needle Insertion

Automatic

Activation Method

Needle Cover

Injection Feedback

Audible click at beginning & end of injection, viewing window

Handling Instructions



1
Remove cap



2
Press & hold against injection site to activate



+Discard in sharps container

Key Features

- ✓ Simple 2-step operation
- ✓ Needle cover activation
- ✓ Automatic needle insertion
- ✓ Audible & visual feedback
- ✓ One-handed operation
- ✓ Ultra-compact design
- ✓ Permanently hidden needle
- ✓ Passive locking needle cover
- ✓ Preconfigured or customized to client specifications
- ✓ Compatible with Rigid Needle Shield