



Maggie[®]

Two Step Cartridge-Based Solution

Volume 

Viscosity 

Device Type

Disposable Auto Injector

Primary Container

3mL Standard ISO Cartridge with Pre-Installed Needle (NIT[®])

Dosage Type

Fixed Dose

Route of Administration

Subcutaneous | Intramuscular

Needle Insertion

Manual

Activation Method

Needle Cover

Injection Feedback

Continuous clicks during injection, viewing window

Handling Instructions



1 Hold device with the cap upright, twist & remove cap



2 Press & hold against injection site to activate



+ Discard in sharps container

Key Features

- ✓ Simple 2-step operation
- ✓ Needle cover activation
- ✓ Manual needle insertion
- ✓ Continuous audible clicks & visual feedback
- ✓ One-handed operation
- ✓ Permanently hidden needle
- ✓ Compact design
- ✓ Passive locking needle cover
- ✓ Customized to client specifications