



Maggie®

Two Step Cartridge-Based Solution

Device Type	Disposable Auto Injector
Primary Container	3mL Standard ISO Cartridge with Pre-Installed Needle (NIT®)
Dosage Type	Fixed Dose
Route of Administration	Subcutaneous Intramuscular
Needle Insertion	Manual
Activation Method	Needle Cover
Injection Feedback	Continuous clicks during injection, viewing window

Handling Instructions



1 Hold device with the cap upright, twist & remove cap



2 Press & hold against injection site to activate



3 Discard in sharps container

Key Features

- ✓ Simple 2-step operation
- ✓ Needle cover activation
- ✓ Manual needle insertion
- ✓ Continuous audible clicks & visual feedback
- ✓ One-handed operation
- ✓ Permanently hidden needle
- ✓ Compact design
- ✓ Passive locking needle cover
- ✓ Customized to client specifications