



**Molly**®  
1.0 mL  
The ideal safety solution

Volume 

Viscosity 

Device Type	Disposable Auto Injector
Primary Container	1mL Long Syringe with Staked Needle   FNS
Dosage Type	Fixed   Full Contents
Route of Administration	Subcutaneous
Needle Insertion	Manual
Activation Method	Needle Cover
Injection Feedback	Audible click at beginning & end of injection, viewing window

## Handling Instructions



Remove cap



Press & hold against injection site to activate



Discard in sharps container

## Key Features

- ✓ Simple 2-step operation
- ✓ Needle cover activation
- ✓ Manual needle insertion
- ✓ Audible, visual & tactile feedback
- ✓ One-handed operation
- ✓ Ultra-compact design
- ✓ Permanently hidden needle
- ✓ Passive locking needle cover
- ✓ Preconfigured or customized to client specifications