



**Molly<sup>®</sup> RNS**  
The ideal safety solution

Volume 

Viscosity 

<b>Device Type</b>	Disposable Auto Injector
<b>Primary Container</b>	1mL Long Syringe with Staked Needle
<b>Dosage Type</b>	Fixed   Full Contents
<b>Route of Administration</b>	Subcutaneous
<b>Needle Insertion</b>	Manual
<b>Activation Method</b>	Needle Cover
<b>Injection Feedback</b>	Audible click at beginning & end of injection, viewing window

## Handling Instructions



Remove cap



Press & hold against injection site to activate



Discard in sharps container

---

## Key Features

- ✓ Simple 2-step operation
- ✓ Needle cover activation
- ✓ Manual needle insertion
- ✓ Audible, visual & tactile feedback
- ✓ One-handed operation
- ✓ Ultra-compact design
- ✓ Permanently hidden needle
- ✓ Passive locking needle cover
- ✓ Preconfigured or customized to client specifications