



Molly®
The ideal safety solution

Volume 

Viscosity 

Device Type	Disposable Auto Injector
Primary Container	2.25mL Long Pre-Filled Syringe RNS Compatible
Dosage Type	Fixed Full Contents
Route of Administration	Subcutaneous
Needle Insertion	Manual
Activation Method	Needle Cover
Injection Feedback	Audible click at beginning & end of injection, viewing window

Handling Instructions



Remove cap



Press & hold against injection site to activate



Discard in sharps container

Key Features

- ✓ Simple 2-step operation
- ✓ Needle cover activation
- ✓ Manual needle insertion
- ✓ Audible, visual & tactile feedback
- ✓ One-handed operation
- ✓ Larger volume
- ✓ Permanently hidden needle
- ✓ Ergonomic pull- or twist-off cap
- ✓ Robust anti-roll design