



SDI MIX[®]+NIT[®]

Twist & Mix

Volume 

Viscosity 

Device Type	Disposable Auto Injector
Primary Container	1mL Dual Chamber Cartridge with Pre-Installed Needle (NIT [®])
Dosage Type	Fixed Full Contents
Route of Administration	Subcutaneous
Needle Insertion	Manual
Activation Method	Needle Cover
Injection Feedback	Audible click at beginning of injection, viewing window

Handling Instructions



1 Turn rear dial to mix



2 Hold device with the cap upright. Twist & remove cap. Needle cover will extend & device will automatically prime



3 Press & hold against injection site to activate



+ Discard in sharps container

Key Features

- ✓ Simple 3-step operation
- ✓ Needle cover activation
- ✓ Manual needle insertion
- ✓ Audible, visual & tactile feedback
- ✓ Needle Isolation Technology (NIT®)
- ✓ Automatic reconstitution & priming
- ✓ Permanently hidden needle